MGUARDPrime
Embolic Protection Stent

Instructions for Use

Table of contents

No. Language

1.	English.....	3
2.	German.....	7
3.	Spanish.....	12
4.	French.....	16
5.	Italian.....	20
6.	Polish.....	24
7.	Hungarian.....	29
8.	Russian.....	33
9.	Turkish.....	38
10.	Greek.....	42
11.	Dutch.....	47
12.	Swedish.....	51
13.	Finnish.....	55
14.	Danish.....	59
15.	Portuguese.....	63
16.	Latvian.....	67
17.	Slovenian.....	71

Instructions for Use - English

1. Description

MGuard Prime™ coronary Stent System is a coronary balloon expandable stent, wrapped with a fiber sleeve, and premounted on a balloon catheter

MGuard Prime™ stent platform is a laser cut cobalt chromium alloy (CoCr) bare stent. It is constructed of circumferential strut rings that are connected to each other with flexible links.

The MGuard Prime™ sleeve is knitted from a Polyethylene Terephthalate (PET) microfiber that creates a flexible mesh cells with a diameter of 180 micron.

The Stent Delivery System Catheter is a double lumen catheter for rapid exchange (RX) with a semi-compliant balloon. Two radiopaque markers, proximal and distal to aid the balloon positioning under fluoroscopy. There is a Luer lock connector (hub) for balloon inflation. MGuard Prime™ is supplied in specific diameters and lengths (see chart and labeling). Nominal stent length and diameter are also printed on the catheter hub.

2. Indications

MGuard Prime™ Coronary Stent System is indicated for improving luminal diameter in vessels with reference diameter from 2.5 to 4.0 mm having lesion length <38 mm and providing embolic protection in the following cases:

- Symptomatic coronary artery disease due to discrete de novo or restenotic lesion in native coronary artery.
- Symptomatic coronary artery disease due to culprit lesion in saphenous vein graft.
- Treatment of coronary lesion in patients undergoing primary or rescue PCI for acute ST-segment elevation myocardial infarction (STEMI)
- Treatment of coronary lesion having atherothrombotic appearance in patients with non ST-elevation acute coronary syndromes (unstable angina and non ST-segment elevation myocardial infarction (NSTEMI)).

3. Contraindications

General contraindications for coronary stenting and for the use of this device are:

- Unprotected coronary left main artery disease.
- Coronary stenosis on a bifurcation.
- Any patients judged to have a lesion which may prevent correct stent deployment.
- Contraindication for anti-platelet and/or anti-coagulation therapy.
- Known allergies to contrast media, Co-Cr or Polyethylene Terephthalate.

4. Warnings

- PTCA and coronary stent implantation are recommended at hospitals where emergency bypass graft surgery can be performed in the event of a potential injury or life-threatening complication. A surgical team should be on standby when an interventional procedure is being performed.
- Only physicians thoroughly trained in percutaneous trans luminal coronary angioplasty (PTCA) and coronary stenting should use this device.
- Cautious selection of patients is necessary since the use of this device carries the associated risk of sub-acute thrombosis, vascular complications and/or bleeding events.
- This device is designed and intended for single use only. DO NOT REPROCESSE, RESTERILIZE AND/OR REUSE. InspireMD will not be responsible for any direct, incidental or consequential damages resulting from reprocessing, re-sterilization or reuse. DO NOT re-sterilize and/or reuse it, as this can potentially result in compromised device performance and risk of cross contamination.
- Check expiry date before using the device. If expiry date has passed, do not use the device.
- Never implant stents of different material close together or

Stent Diameter (mm)	Stent Lengths Available (mm)	Minimum Guiding catheter compatibility	Maximum guidewire compatibility	Nominal stent deployment pressure (atm)	Rated burst pressure (atm)	Strut thickness (micron)
2.5	8,13,18,23,28,33,38	6F (1,78mm)	0,014"(0,36mm)	8	16	80
2.75	8,13,18,23,28,33,38	6F (1,78mm)	0,014"(0,36mm)	8	16	80
3.0	8,13,18,23,28,33,38	6F (1,78mm)	0,014"(0,36mm)	8	16	80
3.25	8,13,18,23,28,33,38	6F (1,78mm)	0,014"(0,36mm)	8	16	80
3.5	8,13,18,23,28,33	6F (1,78mm)	0,014"(0,36mm)	8	16	80
3.5	38	6F (1,78mm)	0,014"(0,36mm)	8	14	80
4.0	8,13,18,23,28	6F (1,78mm)	0,014"(0,36mm)	8	16	80
4.0	33,38	6F (1,78mm)	0,014"(0,36mm)	8	14	80

overlapping the struts (due to galvanic potential risk).

- Should any resistance be felt at any time during manipulation, insertion or withdrawal of the device, do not force: stop the procedure immediately and determine the cause of the resistance before proceeding. If the stent cannot be deployed, remove the entire stent system and guiding catheter as a unit, do not attempt to retrieve the stent back through the guiding catheter dislodgement of the stent may occur.
- Appropriate anti-platelet and/or anti-coagulant therapy should be administered to the patient as determined by the physician in accordance with standard protocols for stent implantation.
- When treating multiple lesions, the distal lesion should be stented initially, followed by stenting of the proximal lesion. Stenting in this order eliminates the need to re-cross the proximal stent and reduces the chance for stent dislodgment.
- Placement of the stent has the potential to compromise side branch patency. Side branches with diameter of 2.0 mm and more should not be covered by the stent.
- Safety and effectiveness of MGuard Prime™ has not yet been established in pregnant patients or patients under the age of 18.

5. Precautions

- Select correct device size suiting lesion length and reference vessel diameter.
- Do not attempt to reposition a partially deployed stent as this may result in severe vessel damage.
- Do not attempt to remove, readjust or touch the stent on its delivery system. The stent cannot be removed and placed on another balloon catheter.
- Balloon pressure should not exceed the Rated Burst Pressure (RBP). Use of a pressure-monitoring device is mandatory to prevent over-pressurization.
- Use only an appropriate balloon inflation medium (a 50:50 mixture by volume of contrast medium and sterile saline). Never use air or a gaseous medium to inflate the balloon.
- When the catheter is in the body, it should be manipulated while under high quality Fluoroscopy.
- While removing the stent protection cap, hold the cap at the most distal part. Holding the cap in the proximal part may cause stent dislodgment.
- Use only guidewires with a maximum of 0.014" (0.36 mm) diameter.
- Exercise care during handling in order to avoid possible damage to the stent, sleeve, balloon-catheter and stent-balloon fixation. Avoid acute bending or kinking of the catheter shaft.
- Store in a dark, cool and dry place.
- Inspect the device prior to the procedure, to verify functionality and lack of damaged parts. Do not use the device if the outer or the inner package is damaged or opened.

6. Complications

Possible complications associated with coronary stent implantation may include, but are not limited to:

- Death
- Acute myocardial infarction
- Coronary artery spasm
- Restenosis of the stented artery
- Acute or sub-acute stent thrombosis
- Total occlusion of the coronary artery or bypass graft
- Dissection, perforation or rupture of the coronary artery
- Hypo/hypertension
- Sepsis/infection
- Arrhythmia, including ventricular fibrillation
- Hemorrhage or hematoma
- Arteriovenous fistula
- Embolization
- Stent migration
- Allergic reactions
- Stent Collapse
- Stent breakage or fracture may occur during implantation
- Drug reactions to anti-platelet agents / contrast medium
- Emergent Coronary Artery Bypass Surgery
- Pain at insertion site
- Myocardial Ischemia,
- Femoral Pseudo-aneurysm,
- Stent thrombosis
- Stroke / cerebrovascular accident
- Renal failure

7. Directions for use

Caution: If any resistance is felt during lesion crossing or if the system is unable to cross the lesion easily, do not attempt to pull the stent system back through the guiding catheter: dislodgement of the stent can occur. Pull back carefully the entire system (guiding catheter, guidewire and stent delivery system) as a single unit. The lesion then must be predilated or otherwise prepared via an alternative strategy, if intention for stenting is maintained. Do not use the same stent system after predilatation because the stent and/or the stent delivery system may have been damaged during the first attempt to cross the lesion or during withdrawal.

Caution: Special care must be taken not to handle or in any way disrupt the stent position on the delivery system. This is most important during catheter removal from packaging, removing the protective cap, placement over guidewire, and advancement through rotating haemostatic valve adapter and guiding catheter hub. Manipulation, e.g., rolling the mounted stent, may cause dislodgement of the stent from the delivery balloon.

7.1. Stent Delivery System Preparation

1. using sterile technique remove the Stent Delivery System and dispenser from the package and place it onto the sterile field.
2. Gently pull out the catheter from the dispenser.
3. Carefully remove the protective cap.
4. Visually check the stent crimping for uniformity, no protruding struts and centering on the balloon.
5. Visually check the sleeve for uniformity, no run outs and centering on the stent. Do not use the device if any defect is noted.

7.2. Flushing of Guidewire Lumen

1. Using a syringe containing sterile saline solution to the distal tip of the catheter, flush the guidewire lumen.
2. Remove the syringe.
3. Avoid manipulation of the stent during removal from packaging and flushing of guidewire lumen.
4. Do not pre-inflate the balloon prior to stent deployment. Use balloon purging technique described further.

Caution: Do not apply negative pressure to the catheter prior to placement of the stent across the lesion. This may cause premature dislodgement of the stent.

7.3. Insertion

1. Visually inspect the stent placement prior to insertion through the hemostatic valve.
2. Attach haemostatic valve to the Luer lock connector of the guiding catheter positioned within the vasculature.
3. In case the coronary guidewire has not already been positioned across the lesion, under fluoroscopy insert a 0.014" (0.36 mm) guidewire across the lesion, following standard PTCA techniques.
4. Insert proximal end of the guidewire into the distal tip of the catheter.
5. Carefully insert the Stent Delivery System through the haemostatic valve and advance the catheter.

Caution: Make sure that the haemostatic valve is completely open before inserting the delivery system, since a partially open valve might damage the stent or dislodge it from the centered location on the dilation balloon.

6. Under fluoroscopic guidance, carefully advance the stent and the Stent Delivery System through the guiding catheter and into the coronary artery following the guidewire towards the target lesion.
7. Position the stent within the lesion using the radiopaque markers located on the balloon catheter as reference points for precise placement across the target lesion.

Caution: If any resistance is felt when the catheter is pushed out of the end of the guiding catheter, be sure to determine the cause. If you cannot determine and eliminate the cause, then retract the entire system (guidewire, delivery system and guiding catheter) together, as a single unit.

7.4. Stent deployment and Balloon removal

1. Maintaining the stent within the lesion, connect a syringe containing 50:50 mixture by volume of contrast medium and sterile saline solution to the Luer lock in the proximal hub of the catheter and apply negative pressure for about 15 seconds until no bubbles appear in the contrast medium solution. Return the syringe leaving a meniscus of contrast in the hub of the balloon lumen.
2. Remove air bubbles from the inflation device following the manufacturer's instructions.
3. Attach inflation device to the SDS. Avoid allowing air to enter the system.
4. Inflate the dilation balloon gradually to expand the stent to the calculated diameter in accordance with the Compliance Chart. Apply constant pressure for at least 30 seconds.
5. Stent should be expanded to a vessel/stent diameter ratio of no more than 1:1.15.

Caution

- Do not expand the stent if it is not properly positioned in the target lesion.
- Stent placement may compromise side branch patency.
- Do not exceed Rated Burst Pressure (RBP).
- In case stent size is still inadequate (under-dilated) with respect to vessel diameter, in order to achieve optimal stent deployment inflate the balloon once more with the same delivery system or perform a post dilatation using a non-compliant high-pressure balloon. Use multiple fluoroscopy views to ensure that the stent has been completely expanded.

7.5. Balloon deflation and catheter removal

1. Deflate the balloon in accordance with standard PCTA procedures. Apply negative pressure to the balloon for about 30 sec. Ensure balloon is fully deflated before pulling the catheter carefully out of the vessel.
2. Observation of the patient and angiographic evolution should be performed periodically in the 15 minutes after stent implantation.

8. MRI Safety Information

The MGuard Prime™ stent was determined to be MRI-conditional through non-clinical testing. A patient with this implant can be scanned safely, immediately after placement under the following conditions:

1. Static magnetic field of 3 Tesla or less
 2. Maximum spatial gradient magnetic field of 720 Gauss/cm or less
- The maximum temperature rise during MRI that was measured was 2.2°C.

MR image quality may be compromised if the area of interest is in exact same area or relatively close to the position of the MGuard Prime™ stent. Therefore optimization of MR imaging parameters to compensate for the presence of this implant may be necessary

9. Pregnancy

This product has not been tested in pregnant women, women who want to be pregnant, nursing mothers or in men intending to father children, effects on the developing fetus have not been studied. While there is no contraindication, the risks and reproductive effects remain unknown.

10. Recommended storage

Store at room temperature and in a dry and dark place. For exact values see labeling.

11. Disposal instructions

After use, dispose the product and packaging in accordance with hospital, administrative and/or local government policy.

12. Reprocessing

This MGuard Prime™ is designed and intended for single use only. Reprocessing for secondary use may degrade the device integrity, sterility and functionality and put the patient in risk. InspireMD will not be responsible for any direct, incidental or consequential damages resulting from reprocessing, resterilization or reuse.

13. Warranty/Liability

The product and each component of its system have been designed, manufactured, tested and packaged with all reasonable care. The warnings contained in InspireMD instructions for use are expressly considered as an integral part of this provision. InspireMD warranties the product until the expiration date indicated on the same. The warranty is valid provided that the use of the product was consistent with the instructions for use. InspireMD disclaims any warranty of merchantability or fitness for a particular purpose of the product. InspireMD is not liable for any direct, indirect, incidental or consequential damages caused by the product. Except in the case of fraud or grave fault on InspireMD part, compensation of any damage to the buyer will not, in any event, be greater than the invoice price of the disputed products. The guarantee contained in this provision incorporates and substitutes the legal guarantees for defects and

compliance, and excludes any other possible liability of InspireMD, however originating, from its product supplied. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects InspireMD legitimate interest in limiting its liability or warranty. No person has any authority to bind InspireMD to any warranty or liability regarding the product.

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InspireMD Ltd.**